

vrinda-pc:8080/dcm4chee-arc/ui2/#/study/study 14:43:00

Patients Studies Series MWL MPPS Work Items Compare

### Studies

Fuzzy Matching  
 Patient family name Patient ID Issuer of patient  
 Accession number Issuer of accession number  
 Study Description Modality  
 Referring physician family Institutional Department N  
 Sending AET of Series Study Instance UID  
 Study date Study time  
 Study Received Study Access

Order By: DCM4CHEE

SUBMIT COUNT SIZE

More More functions

**Upload DICOM data**

Select the STOW-RS server:  
DCM4CHEE

Choose DICOM Files to upload:

Browse... STROKE^Emergency.dcm

STROKE^Emergency.dcm ✓

CLOSE

vrinda-pc:8080/dcm4chee-arc/ui2/#/study/study 14:43:46

Patients Studies Series MWL MPPS Work Items Compare

### Studies

Fuzzy Matching  
 Patient family name Patient ID Issuer of patient  
 Accession number Issuer of accession number  
 Study Description Modality  
 Referring physician family Institutional Department N  
 Sending AET of Series Study Instance UID  
 Study date Study time  
 Study Received Study Access

Order By: DCM4CHEE

SUBMIT COUNT SIZE

More Search current list More functions Actions for selections

	Patient's Name	Patient ID	Issuer of Patient	Birth Date	Sex	Patient Comments	#S
	Emergency STROKE	J4C-UC-001		19390825	F		1
1	UC-UPS-	1.3.12.2.1107.5.8.3.485251.834954.83838049.2020	20200910 054928 Neuro 8/4	043J4C2020UC001	1620101352	MR HEAD^STANDARD	1 1

vrinda-pc:8080/dcm4chee-arc/ui2/#/device/edit/dcm4chee-arc 14:44:19

device list dcm4chee-arc

Search in this part...

- Extensions
- Child Objects
- Attributes

<b>Device Name</b>	dcm4chee-arc	A unique name for this device
<b>Device Description</b>		Unconstrained text description of the device
<b>Vendor Device Data</b>	Download	Device specific vendor configuration information
<b>Device UID</b>		Unique identifier of the device
<b>Manufacturer</b>	dcm4chee.org	Should be the same as the value of Manufacturer (0008,0070) in SOP instances created by this device
<b>Manufacturer Model Name</b>	dcm4chee-arc	Should be the same as the value of Manufacturer Model Name (0008,1090) in SOP instances created by this device
<b>Software Version</b>	5.24.2 ✕ b8b0c20 ✕ 2021-10-01 ✕ + Add	Should be the same as the values of Software Versions (0018,1020) in SOP instances created by this device